DICOM media startup guide **(ENGLISH)**

This short guide explains how to launch the viewer application which allows displaying the DICOM medical images present on this portable media (CD, DVD, USB stick, etc.). Depending on the configuration of the application which created this media, the viewer application may also be absent from the present media. Instructions for each supported operating system follow.

Windows operating system

Nowadays auto running of removable media is disabled for security reasons on most Windows systems. Hence, in the likely event the *LocalEye Viewer* application does not start automatically upon media insertion, please browse to the root of the media contents, and double-click on the “*ViewMediaWindows.exe*” file: the viewer application will start.

Please refer to the Online Help of the viewer application for further information about its usage.

macOS operating system

If you are using a Mac, please insert the media, wait a few seconds, then double click on the media icon which will appear on your desktop, in order to open the media itself and browse its contents. Double-click on the “*ViewMediaMac*” file: a new *LocalEyeViewer* volume will be automatically mounted, and its contents will be shown in a new window. At this point, double-click on the *Launch LocalEye Viewer* icon in this window: the viewer application will start.

Please refer to the Online Help of the viewer application for further information about its usage.

Linux operating system

If you are using a Linux PC with a graphic windowing environment, please insert the media, wait a few seconds, then double click on the media icon which will appear on your desktop (auto-mounting shall be enabled), in order to open the media itself and browse its contents. Double-click on the “*ViewMediaLinux.sh*” file: the viewer application will start.

Please refer to the Online Help of the viewer application for further information about its usage.

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